

HD IntegraTM

DRS-1 Dialyzer Rinsing System



Features:-

- Color touch screen display
- 4 channel rinsing station
- Quick rinsing process (approx. 5mins)
- Digital monitoring of disinfectant in dialyzer during rinsing process
- 4 channel rinsing status display
- Automatic calibration for RO water conductivity
- Digital water inlet pressure monitoring
- Save on saline usage
- One touch Sanitize mode
- Easy maintenance
- Conforms to IEC60601 safety requirements
- Compact design
- Audio and visual alarm display

Introduction

DRS-1 is developed to save time and saline usage in rinsing of reprocessed dialyzers. The system efficiently removes disinfectant inside the dialyzer without contaminating the blood compartment.

The system utilizes digital sensors to monitor dialyzer disinfectant level to ensure disinfectant is completely removed before use. This feature adds safety in the rinsing process, reduces water usage and reduces rinsing time.

To ensure safety, the system provides automatic Sanitization to the system. Any commercial disinfectant e.g. Bleach and PAA is suitable for sanitizing the system.

FSA-TECHNOLOGY INC.

Unit 2ARLI Bldg., South Point, Banay-banay
Cabuyao City Laguna, Philippines 4025
Tel #: +63 49 5446972 Email: fsasales@fsa-tech.com